

Dust test equipment

Description:

The dust chamber, as specified in the standard, incorporates a blower, circulating ducting, and a standard headform.

A suitable reflectometer head and digital display of reflectance are included.

An adequate supply of the required dust is also supplied.

Services required:

Floor mounted
110/230 volts AC, 50/60Hz, mains electricity

Approximate packed size & weight:

200 x 60 x 80 cm : 100 kg

Relevant standards:

EN 168:2001, clause 13.1 and figure 6
ANSI Z87.1:2010, clause 9.17.2
AS/NZS 1337-1:2010, appendix W

