

Gas test equipment

Description:

The equipment comprises a steel cabinet, as specified, fitted with an inlet and a fan assisted outlet for the Ammonia gas.

The test headform, mounted on a wooden base, is included.

The required gas washing bottle, hand bellows and connecting tubing (not shown in the photo) are also supplied.

Services required:

Bench mounted
110/230 volts AC, 50/60Hz, mains electricity
Exhaust of Ammonia gas to outside
Supply of concentrated liquid Ammonia
Supply of detecting solution chemicals specified in standard

Approximate packed size & weight:

80 x 60 x 60 cm, 25 kg

Relevant standards:

EN 168:2001, clause 14.1
ANSI Z87.1:2010, clause 9.18.2
AS/NZS 1337-1:2010, appendix X

