

## European reference headforms (medium and small)

### **Description:**

These headforms are made of cast epoxy resin (finished white painted) and engraved.

Each headform is provided with two removable 'ear pins' to support spectacle side arms.

These headforms are used as reference heads for evaluating eye and face protectors for areas to be protected, peripheral vision, impact points, etc.

### **Services required:**

None

### **Approximate packed size & weight:**

40 x 35 x 25 cm : 7 kg

### **Relevant standards:**

EN 168:2001, clause 17 and figure 11  
AS/NZS 1337-1:2010, clause 1.4.6 and figure 1

*Medium headform*



*Small headform*

