

Illuminated grid test equipment

Description:

An illuminated grid on a diffuse viewing area, approximately 30cm x 20cm.

The grid is viewed from a distance of approximately 1m, looking through the ocular being assessed.

Any irregularities or gross imperfections in the ocular are easily observed.



Services required:

Bench or wall mounted
110/230 volts AC, 50/60Hz, mains electricity

Approximate packed size & weight:

50 x 40 x 20 cm, 3 kg

Relevant standards:

EN 167:2001, clause 5
ANSI Z87.1:2010, clause 9.1.2
AS/NZS 1337-1:2010, appendix J