

Acoustic test fixture (ATF)

Description:

As fully described in the standards, and required for testing earmuffs to EN 352 Parts 1 and 3, the acoustic test fixture is designed to accommodate a 1" or 1/2" diameter microphone and pre-amplifier.

It is provided with an isolation cup and (optionally) an adaptor to accommodate earmuffs mounted on an industrial safety helmet.

Services required:

None

Approximate packed size & weight:

40 x 30 x 30 cm : 8 kg

Relevant standards:

EN 13819-2:2002, clause 4.1.2
ISO 4869-3:2007, clause 5.1 and figure 1

