

Drop test plate

Description:

As required for testing earmuffs to in EN 352 Parts 1 and 3, and earplugs to Part 2, this mild steel plate, 500 x 500 x 12 mm, is edged & chamfered, flash ground top and bottom and zinc plated.

To be placed on a firm, level floor.

Services required:

Floor mounted

Approximate packed size & weight:

70 x 70 x 30 cm : 7 kg

Relevant standards:

EN 352 – required for testing to all parts of this standard
EN 13819-1:2002, clauses 4.6.2, 5.4.2 and 5.5.2.2

