

Temperature controlled water bath

Description:

As specified in the standard and required for conditioning earmuffs to in EN 352 Parts 1 and 3, this stainless steel bath is designed to accommodate several specimens (headband or helmet mounted earmuffs).

The temperature controller maintains the water temperature within the specified range of 50 +/-2°C.

Services required:

Bench mounted
110/230 volts AC, 50/60Hz, mains electricity

Approximate packed size & weight:

100 x 60 x 50 cm : 25 kg

Relevant standard:

EN 13819-1:2002, clause 4.10

