

Stolle apertometer

Description:

The apertometer comprises two major assemblies:

A headform with light bulbs located at the eye positions and;

A graduated, translucent, thermo-formed plastic bowl of nominal 610mm diameter.

These are fixed to a suitable frame to enable bench mounting.

Also included is a planimeter, to measure the areas of the field of vision chart, produced using the apertometer.



Services required:

Bench mounted
110/230 volts AC, 50/60Hz, mains electricity

Approximate packed size & weight:

90 x 60 x 50 cm : 25 kg

Relevant standard(s):

EN 136:1998, clause 8.17 and figure 14
EN 250: 2006, clause 6.8.2.3 and figure 2
EN 12941, clause 7.4 and figure 3
EN 12942, clause 6.3.2.4.2