

Filter Simulators

Description:

The filter simulator is a device that has a specified mass, size and pressure drop according to the relevant standard.

There are three versions which cover European and Australian standards. They each incorporate an EN 148-1 male thread.



Available versions:

500g, 96mm (AS/NZS 1716:2012, EN 136:1998)

300g, 96mm (EN 140:1998)

300g, 50mm (EN 142:2002)

Services required:

Breathing machine or breathing resistance test equipment or inward leakage test equipment as required

Approximate packed size & weight:

45 x 10 x 10cm (set of three), 1.5kg

Relevant standards:

AS/NZS 1716:2012

EN 136:1998

EN 140:1998

EN 142:2002

EN 13274-1:2001