

Supplementary CO₂ test equipment (for use with breathing machine)

Description:

The equipment set comprises a single channel CO₂ analyser, flow meter, non-return valves, compensators and a CO₂ absorber.

A second CO₂ analyser can be supplied to provide simultaneous measurements of both the 5% CO₂ concentration and sample CO₂ concentration.



Services required:

Bench mounted
110/230 volts AC, 50/60Hz, mains electricity
CO₂ supply cylinder and regulator
Reference CO₂ gases
Additional air-flow duct
Sheffield head/torso
Breathing machine

Approximate packed size & weight:

To be advised with quotation

Relevant standards:

EN 136
EN 140, EN 149
EN 12941, EN 12942
EN 13274-3, EN 13274-6
EN 14594