

## Polarizing axis test equipment

### Description

The equipment comprises a pair of individually mounted split field polarizers, joined and glass mounted. These can then be rotated, the degree of rotation being displayed by a pointer passing over a calibrated scale.

The polarizers are illuminated from behind by a diffuse light source.

The method of use is given in the standards.



### Services required

Bench mounted  
110/230 volts AC, 50/60Hz, mains electricity

### Approximate size & weight

270 x 170 x 120 mm : 1 kg

### Relevant standards

EN 172:2001, annex D  
EN 1836:2007, clause 6.2.3.3.2  
AS/NZS 1067.1:2016, clause 5.3.4.2  
AS/NZS 1067.2:2016, clause 7.10.1  
AS/NZS 1337-1:2010, appendix C  
ISO 12311:2013, clause 7.10.1  
ISO 12312:2013, Part 1, clause 5.3.4.2