

## Set of test headforms for helmet-mounted earmuffs

### **Description:**

Each of the three headforms (size codes B, J and N) comprises a cast epoxy resin “half-head”, secured between two parallel, transparent acrylic plates.

On each plate, at the points corresponding to the centre of the ear openings, is located a threaded insert to enable the attachment of an artificial pinna at the appropriate point.



The headform and plates assembly is attached to a hardwood baseplate.

One pair of the specified artificial pinnae and securing screws are included.

### **Services required:**

Bench mounted

### **Approximate packed size & weight:**

80 x 40 x 35 cm : 10 kg

### **Relevant standard:**

EN 13819-1:2002, clauses 4.2.2.2 and figure 4