

Stolle apertometer

Description

The apertometer comprises two major assemblies:

A headform with light bulbs located at the eye positions and;

A graduated, translucent, thermo-formed plastic bowl of nominal 610mm diameter.

These are fixed to a suitable frame to enable bench mounting.



Services required

Bench mounted
110/230 volts AC, 50/60Hz, mains electricity

Optional items

A planimeter (e.g. 9PL020S), or similar, is required to measure the areas of the field of vision chart, produced using the apertometer.

Approximate sizes & weights

Apertometer Bowl: 660 x 520 x 690 mm : 10kg
Apertometer Head: 170 x 210 x 460 mm : 5kg

Relevant standards

EN 136:1998, clause 8.17 and figure 14
EN 250:2014, clause 6.10.2.3 and figure 2
EN 12941, clause 7.4 and figure 3
EN 12942, clause 6.3.2.4.2

Please [contact us](#) for a quotation

1AP1909