

## Filter test chamber

### Description

The filter test chamber is an add-on to a NaCl flame photometer, allowing filters and filtering facepiece to be tested. It can also be used for filter gas testing.

The test chamber consists of a clear body with aluminium bulkheads. A locking clamp ensures that the door closes tightly against the door seal, preventing any leakage.



### Services required

Bench mounted.

### Approximate size & weight

50x 40 x 40 cm : 10 kg

### Relevant standards

- AS 1716
- EN 143
- EN 149
- EN 403
- EN 405
- EN 1827
- EN 12941
- EN 12942
- EN 13274-7
- EN 14387

Please [contact us](#) for a quotation

Issue: 1FTC.18.07

Your global partner for testing and certification of all types of personal protective equipment